

HDGFRP3 Rabbit pAb

货号: B15880

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	22kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse kidney,Mouse lung,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 90-203 of human HDGF RP3 (NP_057157.1).
序列	NNPGVKFTGYQAIQQSSSETEGEGGNTADASSEEEGDRVEEDGKGKRKNEKAGSKRKKSYSKSSKQSRKSPGDED DKDCKEEENKSSSEGGDAGNDTRNTTSDLQKTSEGT

靶点信息

研究背景	Predicted to enable double-stranded DNA binding activity; microtubule binding activity; and transcription coregulator activity. Predicted to be involved in several processes, including microtubule polymerization; negative regulation of microtubule depolymerization; and neuron projection development. Located in nucl eoplasm.
基因ID	50810

基因名	HDGFL3
Swiss	
别名	HDGFRP3;CGI-142;HDGF-2;HDGF2;HRP-3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn