

PBXIP1 Rabbit pAb

货号: B15887

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:1000 - 1:2000 IHC: 1:50 - 1:200
理论分子量	72kDa/77kDa/80kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,LO2,BT-474,A-549
细胞定位	Cytoplasm,Nucleus,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human PBXIP1 (NP_065385.2).
序列	MASCPDSDNSWVLAGSESLPVETLGPASRMDPESERALQAPHSPKTDGKELAGTMDGEGTLFQTESQSGSILTEETE VKGTLEGDVCGVEPPGPGDTVVQGDLDQETTVVTGLGPDQTDLGQSPQSLPSTPKAAWIREEGRCSSDDDDTDVDME GLRRRRRGREAGPPQPMVPLAVENQAGGEGAGGE

靶点信息

研究背景	The protein encoded by this gene interacts with the PBX1 homeodomain protein, inhibiting its transcriptional activation potential by preventing its binding to DNA. The encoded protein, which is primarily cytosolic but can shuttle to the nucleus, also can interact with estrogen receptors alpha and beta and promote the proliferation of breast cancer, brain tumors, and lung cancer. Several transcript variants encoding different isoforms have been found for this gene. More variants exist, but their full-length natures have yet to be determined.
------	--

基因ID	57326
基因名	PBXIP1
Swiss	Q96AQ6
别名	PBXIP1;HPIP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn