

MRFAP1 Rabbit pAb

货号: B15894

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:1000
理论分子量	14kDa
实测分子量	14kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,A-431,SKOV3,MCF7
细胞定位	Cytoplasm,Nucleus,perinuclear region
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human MRFAP 1 (NP_150638.1).
序列	MRPLDIVELAEPEEEVLEPEEDFEQFLLPVINEMREDIASLTREHGRAYLRNRSKLWEMDNMLIQIKTQVEASEESALNHLQ NPGDAAEGRAAKRCEKAEEKAKEIAKMAEMLVELVRRIEKSESS

靶点信息

研究背景	This gene encodes an intracellular protein that interacts with members of the MORF4/MRG (mortality fact or on chromosome 4/MORF4 related gene) family and the tumor suppressor Rb (retinoblastoma protein.) The protein may play a role in senescence, cell growth and immortalization. Alternative splicing results in multiple transcript variants.
基因ID	93621

基因名	MRFAP1
Swiss	Q9Y605
别名	MRFAP1;PAM14;PGR1

产品验证



Western blot analysis of MRFAP1 expressed in Jurkat,A-431,SKOV3 using MRFAP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn