

SYMPK Rabbit pAb

货号: B15899

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	74kDa/141kDa
实测分子量	141kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,293T,HL-60,HT-1080,A-549,Mouse brain,Mouse liver
细胞定位	Cell junction,Cell membrane,Cytoplasm,Cytoplasmic side,Nucleus,Peripheral membrane protein,cytoskeleton,nucleoplasm,tight junction
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 310-530 of human SYMPK (NP_004810.2).
序列	KHPASLEFQAQITLLVDLGTPQAEIARNMPSSKDTRKRPRDDSDTLKKMKLEPNLGEDDEDKDLEPGPSGTSKASAQIS GQSDTDITAEFLQPLLTPDNVANLVLISMVYLPEAMPASFQAIYTPVESAGTEAQIKHLARLMATQMTAAGLPGPGVEQTKQ CKEKPKEEKVVKTESVLIKRRLSAQGQQAISVVGSLSSMSPLIEEEAPQAKRRPEPIIPVT

靶点信息

研究背景	This gene encodes a nuclear protein that functions in the regulation of polyadenylation and promotes gene expression. The protein forms a high-molecular weight complex with components of the polyadenylation machinery. It is thought to serve as a scaffold for recruiting regulatory factors to the polyadenylation complex. It also participates in 3'-end maturation of histone mRNAs, which do not undergo polyadenylation. The protein also localizes to the cytoplasmic plaques of tight junctions in some cell types.
------	---

基因ID	8189
基因名	SYMPK
Swiss	Q92797
别名	SYMPK;SPK;SYM;symplekin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn