

MIB1/DIP1 Rabbit pAb

货号: B15925

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	110kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-549,HepG2,Mouse stomach,Mouse lung,Mouse kidney ,Rat testis
细胞定位	Cell membrane,Cytoplasm,centriolar satellite,centrosome,cytoskeleton,microtubule organizing center
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 560-860 of human MIB1 /DIP1 (NP_065825.1).
序列	DSEGDTPLHDAISKRRDDILAVLLEAGADVTITNNNGFNALHHAALRGNPSSAMRVLLSKLPRPWIVDEKKDDGYTALHLAA LNNHVEVAELLVHQGNANLDIQNVNQQTALHLAVERQHTQIVRLLVRAGAKLDIQDKDGDTPLHEALRHHTLSQLRQLQD MQDVGKVDAAWEPKNTLIMGLGTQGAEEKSAASIACFLAANGADLSIRNKKGQSPDLCPDPNLCKALAKCHKEKVSG QVGSRSPPMISNDSETLEECMVCSDMKRDTLFGPCGHIATCSLCSPRVKKCLICKEQVQSR

靶点信息

研究背景	This gene encodes a protein containing multiple ankyrin repeats and RING finger domains that functions as an E3 ubiquitin ligase. The encoded protein positively regulates Notch signaling by ubiquitinating the Notch receptors, thereby facilitating their endocytosis. This protein may also promote the ubiquitination and degradation of death-associated protein kinase 1 (DAPK1).
------	--

基因ID	57534
基因名	MIB1
Swiss	Q86YT6
别名	MIB1;DIP-1;DIP1;LVNC7;MIB;ZZANK2;ZZZ6

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn