

HERC3 Rabbit pAb

货号: B15934

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	39kDa/117kDa
实测分子量	117kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,U-87MG,A-431,Mouse lung,Mouse spinal cord
细胞定位	Cytoplasm,Cytoplasmic vesicle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HERC3 (NP_055421.1).
序列	MLCWGYWSLGQPGISTNLQGIVAEPQVCGFISDRSVKEVACGGNHSVFLLLEDGEVYTCGLNTKGQLGHEREGNKPEQIG ALADQHIIHVACGESHSLSLSDRGQLFSWGAGSDGQLGLMTTEDSVAVPRLIQKLNQQTILQVSCGNWHCLALAADGQF FTWGKNSHGQLGLGKEFPSQASPQRVRSLEGIPLAQVAAGGA

靶点信息

研究背景	This gene encodes a member the HERC ubiquitin ligase family. The encoded protein is located in the cyto sol and binds ubiquitin via a HECT domain. Mutations in this gene have been associated with colorectal a nd gastric carcinomas. Alternatively spliced transcript variants encoding multiple isoforms have been obs erved for this gene.
基因ID	8916

基因名	HERC3
Swiss	Q15034
别名	HERC3

产品验证



Western blot analysis of HERC3 expressed in Mouse lung,SH-SY5Y,U-87MG using HERC3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn