

MYO5B Rabbit pAb

货号: B15937

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	47kDa/111kDa/213kDa
实测分子量	213kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,SW480
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1440-1690 of human MYO5B (NP_001073936.1).
序列	AQALAQSERKRHELNRQVTVQRKEKDFQGMLEYHKEDEALLIRNLVTDLKPQMLSGTVPCLPAYILYMCIRHADYTNDLKVHSLLTSTINGIKKVLKKHNDDFEMTSFWLSNTCRLHCLKQYSGDEGFMQTQNTAKQNEHCLKNFDLTEYRQVLSDSLQIYQQLIKIAEGVLQPMIVSAMLENESIQGLSGVKPTGYRKRSSSMADGDNSYCLEAIIRQMNAFHTVMCDQGLDPEIILQVFKQLFYMIN

靶点信息

研究背景	The protein encoded by this gene, together with other proteins, may be involved in plasma membrane re cycling. Mutations in this gene are associated with microvillous inclusion disease.
基因ID	4645

基因名	MYO5B
Swiss	Q9ULV0
别名	MYO5B;myosin VB

产品验证



Western blot analysis of MYO5B expressed in SKOV3,SW480 using MYO5B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn