

C19orf40 Rabbit pAb

货号: B15952

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	23kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human C19orf40 (NP_689479.1).
序列	MEKNPPDDTGPHVPLGHIVANEKWRGSQLAQEMQGKIKLIFEDGLTPDFYLSNRCCILYVTEADLVAGNGYRKRLVRVRNSNNLKGIVVVEKTRMSEQYFPALQKFTVLDLGMVLLPVASQMEASCLVIQLVQEQTKEPSKNPLLGKKRALLLSEPSLLRTVQQIPGVGKVKAPLLLQKFPSIQQLSNASIGELEQVVVGQAVAQQIHAFFTQPR

靶点信息

研究背景	FAAP24 is a component of the Fanconi anemia (FA) core complex (see MIM 227650), which plays a crucial role in DNA damage response (Ciccia et al., 2007 [PubMed 17289582]).[supplied by OMIM, Mar 2008]
基因ID	91442

基因名	FAAP24
Swiss	Q9BTP7
别名	FAAP24;C19orf40

产品验证



Western blot analysis of C19orf40 expressed in MCF7 using C19orf40 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn