

ZNF408 Rabbit pAb

货号: B15957

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	78kDa
实测分子量	78kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HepG2,HeLa,MCF7,Mouse liver,Mouse spleen,Mouse thymus,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 471-720 of human ZNF408 (NP_079017.1).
序列	PVCGRPLANQGSLRNHMLHTGEKPFLCPHCGRAFRQRGNLRGHLRLHTGERPYRCPHCADAFPQLPELRRHLISHTGEAHLCPVCGKALRDPHTLRAHERLHSGERPFPQCGRAYTLATKLRRHLKSHLEDKPYRCPTCGMGYTLPQSLRRHQLSHRPEAPCSPPSVPSAASEPTVVLLQAEPQLLDTHREEEVSPARDVVEVTISESQEKCFVVPEEPDAAPSLVLIHKDMGLGAWAEVVEVEMGT

靶点信息

研究背景	The protein encoded by this gene contains ten tandem zinc fingers and an N-terminal SET domain, so it is likely a DNA binding protein that interacts with other proteins. In adults, the encoded protein is expressed most highly in retina. Consequently, defects in this gene have been associated with familial exudative vitreoretinopathy (FEVR) and retinitis pigmentosa (RP).
------	--

基因ID	79797
基因名	ZNF408
Swiss	Q9H9D4
别名	ZNF408;EVR6;RP72

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn