

NUDT18 Rabbit pAb

货号: B15959

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	24kDa/35kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7
细胞定位	cytosol
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human NUDT18 (NP_079091.3).
序列	MASEGLAGALASVLAGQGSSVHSCDSAPAGEPPAPVRLRKNVCYVVLAVFLSEQDEVLLIQEAKRECRGSWYLPAGRMEPGETIVEALQREVKEEAGLHCEPETLLSVEERGPSWVRFVFLARPTGGILKTSKEADAESLQAAWYPRTSLPTPLRAHDILHLV ELAAQYRQARHPLILPQELPCDLVCQRLVATFTSAQTVWVLVGTVGMPHLPVTACGLDPMEQRGGMKMAVLRLLQECL TLHHLVVEIKGLLGLQHLGRDHS DGICNLVLTVAFRSPGIQDEPPKVRGENFSWWKVMEEDLQSLLQLRLQGSSVVPV NR

靶点信息

研究背景	The protein encoded by this gene is a member of the Nudix hydrolase family. Nudix hydrolases eliminate potentially toxic nucleotide metabolites from the cell and regulate the concentrations and availability of many different nucleotide substrates, cofactors, and signaling molecules. This protein contains a Nudix hydrolase domain and hydrolyzes oxidized forms of guanosine and deoxyguanosine diphosphates.
------	--

基因ID	79873
基因名	NUDT18
Swiss	Q6ZVK8
别名	NUDT18;MTH3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn