

TGFB1I1 Rabbit pAb

货号: B15982

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	47kDa/49kDa
实测分子量	50KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse lung,Rat lung
细胞定位	Cell junction,Cytoplasm,Nucleus matrix,cytoskeleton,focal adhesion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-220 of human TGFB1I1 (NP_001035919.1).
序列	VASGEQKEDQSEDKKRPSLPSSPSPGLPKASATSATLELDRLMASLSDFRVQNHLPASGPTQPPVVSSTNEGSPSPPEPTGKGS�DTMLGLLQSDLSRRGV

靶点信息

研究背景	This gene encodes a coactivator of the androgen receptor, a transcription factor which is activated by an drogen and has a key role in male sexual differentiation. The encoded protein is thought to regulate andr ogen receptor activity and may have a role to play in the treatment of prostate cancer. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	7041

基因名	TGFB1I1
Swiss	O43294
别名	TGFB1I1;ARA55;HIC-5;HIC5;TSC-5

产品验证



Western blot analysis of TGFB1I1 expressed in HeLa,Mouse lung,Rat lung using TGFB1I1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn