

GHRHR Rabbit pAb

货号: B15999

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	47kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Mouse brain,Rat brain
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human GHRHR (NP_000814.2).
序列	MDRRMWGAHVFCVLSPLPTVLGHMHPECDFITQLREDESACLQAAEEMPNTTLGCPATWDGLLCWPTAGSGEWVTLP CPDFFSHFSSGAVKRDCTITGWSEFPFPYPVACPVPLELLAEEESYFSTVK

靶点信息

研究背景	This gene encodes a receptor for growth hormone-releasing hormone. Binding of this hormone to the receptor leads to synthesis and release of growth hormone. Mutations in this gene have been associated with isolated growth hormone deficiency (IGHD), also known as Dwarfism of Sindh, a disorder characterized by short stature.
基因ID	2692

基因名	GHRHR
Swiss	Q02643
别名	GHRHR;GHRFR;GRFR;IGHD1B

产品验证



Western blot analysis of GHRHR expressed in U-251MG,Mouse brain,Rat brain using GHRHR Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn