

NCOR2 Rabbit pAb

货号: B16009

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:100
理论分子量	81kDa/268kDa/273kDa/274kDa
实测分子量	275kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,LO2
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 670-900 of human NCO R2 (NP_006303.4).
序列	KLKMEKERNARRKKKKAPAAASEEAAFPVVEDEEMEASGVSGNEEEMVEEAEALHASGNEVPRGECSGPATVNNSSDT ESIPSPHTEAAKDTGQNGPKPPATLGADGPPPGPPTPPPEDIPAPTEPTPASEATGAPTPPPAPPSPSAPPPVVPKEEKEEET AAAPPVEEGEEQKPPAAEELAVDTGKAEEPVKSECTEEAEEGPAKGKDAAEAEATAEGALKAEKKEGGS

靶点信息

研究背景	This gene encodes a nuclear receptor co-repressor that mediates transcriptional silencing of certain target genes. The encoded protein is a member of a family of thyroid hormone- and retinoic acid receptor-associated co-repressors. This protein acts as part of a multisubunit complex which includes histone deacetylases to modify chromatin structure that prevents basal transcriptional activity of target genes. Aberrant expression of this gene is associated with certain cancers. Alternate splicing results in multiple transcript variants encoding different isoforms.
基因ID	9612
基因名	NCOR2
Swiss	Q9Y618
别名	NCOR2;CTG26;N-CoR2;SMAP270;SMRT;SMRTE;SMRTE-tau;TNRC14;TRAC;TRAC-1;TRAC1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn