

TMEM189 Rabbit pAb

货号: B16011

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100
理论分子量	30kDa/31kDa
实测分子量	31kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7
细胞定位	Endoplasmic reticulum membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human TME M189 (NP_954580.1).
序列	TLLQDWHVILPRKHHRIHHVSPHETYFCITTGWLNYPLEKIGFWRRLEDLIQGLTGEKPRADDMKWAQKIK

靶点信息

研究背景	Co-transcription of this gene and the neighboring downstream gene (ubiquitin-conjugating enzyme E2 variant 1) generates a rare read-through transcript, which encodes a fusion protein comprised of sequence s haring identity with each individual gene product. The protein encoded by this individual gene lacks a UE V1 domain but includes three transmembrane regions. Alternative splicing results in multiple transcript v ariants.
基因ID	387521

基因名	
Swiss	A5PLL7
别名	TMEM189;KUA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn