

KANSL3 Rabbit pAb

货号: B16057

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	16kDa/58kDa/80kDa/83kDa/93kDa/95kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Jurkat
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human KANSL3 (NP_001108488.1).
序列	MAHRGGERDFQTSARRMGTSLLFQLSVHERELDLVFLDHSYAKPWSAHPDASSARPTRMLFVTPRRQHESTIESDVPIDV ETVTSTPMPLYDNQKARSVMNECERHVIFARTDADAPPPPEDWEEHVNRTGWTMAQNKLFNKILKALQSDRLARLANEG ACNEPVLRRVAVDKCARRVRQALASVSWDTKLIQWLHTTLV

靶点信息

研究背景	As part of the NSL complex it is involved in acetylation of nucleosomal histone H4 on several lysine residues and therefore may be involved in the regulation of transcription.
基因ID	55683

基因名	KANSL3
Swiss	Q9P2N6
别名	KANSL3;KIAA1310;NSL3;Rcd1

产品验证



Western blot analysis of KANSL3 expressed in U-251MG, Jurkat using KANSL3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn