

IL36A Rabbit pAb

货号: B16065

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	17kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,BT-474,SW480,A-431,HeLa,Mouse skin,Mouse heart,Rat liver
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-158 of human IL36A (NP_055255.1).
序列	MEKALKIDTPQQGSIQDINHRVWVLQDQTLIAVPRKDRMSPVTIALISCRHVETLEKDRGNPIYLGLNGLNLCLMCAKVGD QPTLQLKEKDIMDLYNQPEPVKSFLFYHSQSGRNSTFESVAFPGWFIASSEGGCPLILTQELGKANTTDFGLTMLF

靶点信息

研究背景	The protein encoded by this gene is a cytokine that can activate NF-kappa-B and MAPK signaling pathway s to generate an inflammatory response. The encoded protein functions primarily in skin and demonstrat es increased expression in psoriasis. In addition, decreased expression of this gene has been linked to a p oor prognosis in both hepatocellular carcinoma and colorectal cancer patients.
基因ID	27179

基因名	IL36A
Swiss	Q9UHA7
别名	IL36A;FIL1;FIL1(EPSILON);FIL1E;IL-1F6;IL1(EPSILON);IL1F6

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn