

SRM Rabbit pAb

货号: B16075

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	33kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	22Rv1,HepG2,HeLa,Mouse liver,Mouse stomach,Mouse small intestine,Rat liver
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-302 of human SRM (N P_003123.2).
序列	MEPGPDGPAASGPAAIREGWFRETCSLWPGQALSQVEQLLHRRSRYQDILVFRSKTYGNVLVLDGVIQCTERDEFSYQ EMIANLPLCSHIPNPRKVLIIGGGDGGVLREVVKHPSVESVVQCEIDEDVIQVSKKFLPGMAIGYSSSKLTLHVGDGFEFMK QNQDAFDVIIITDSSDPMGPAESLFKESYYQLMKTALKEDGVLCCQGECQWLHLDLIKEMRQFCQSLFPVVAYAYCTIPTYP SGQIGFMLCSKNPSTNFQEPVQPLTQQQVAQMQLKYNSDVHRAAFVLPEFARKALNDVS

靶点信息

研究背景	The polyamines putrescine, spermine, and spermidine are ubiquitous polycationic mediators of cell growth and differentiation. Spermidine synthase is one of four enzymes in the polyamine-biosynthetic pathway and carries out the final step of spermidine biosynthesis. This enzyme catalyzes the conversion of putrescine to spermidine using decarboxylated S-adenosylmethionine as the cofactor.
------	---

基因ID	6723
基因名	SRM
Swiss	P19623
别名	SRM;PAPT;SPDSY;SPS1;SRML1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn