

CD6 Rabbit pAb

货号: B16086

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	58kDa/60kDa/63kDa/64kDa/68kDa/71kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Jurkat,U-937,HT-29,Raji
细胞定位	Cell membrane,Secreted,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 424-668 of human CD6 (NP_006716.3).
序列	LRIKGKYALPVMVNHQHLPTTIPAGSNSYQVPVITIPKEVFMPIQVQAPPPEDSDSGSDSDYEHYDFSAQPPVALTTFYNS QRHRVTDEEVQQSRFQMPPLEEGLEELHASHIPTANPGHCITDPPSLGPQYHPRSNSSESSTSSGEDYCNSPKSKLPPWNP QVFSSERSSFLEQPPNLELAGTQPAFSAGPPADDSSSTSSGEWYQNFQPPPQPPSEEQFGCPGSPSPQPDSTDNDDYDD ISAA

靶点信息

研究背景	This gene encodes a protein found on the outer membrane of T-lymphocytes as well as some other immune cells. The encoded protein contains three scavenger receptor cysteine-rich (SRCR) domains and a binding site for an activated leukocyte cell adhesion molecule. The gene product is important for continuation of T cell activation. This gene may be associated with susceptibility to multiple sclerosis (PMID: 19525953 , 21849685). Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	923
基因名	CD6
Swiss	P30203
别名	CD6;TP120

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn