

ADAMTSL1 Rabbit pAb

货号: B16093

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	27kDa/48kDa/50kDa/58kDa/75kDa/193kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,U-251MG,U-87MG,Mouse lung,Mouse kidney
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1483-1762 of human ADAMTSL1 (NP_001035362.3).
序列	SLVIQDYWWSVDRLATCSASCGNRGVQQPRLRCLLNSTEVNPAHCAGKVRPAVQPIACNRRDCPSRWMVTSWSACTRS CGGGVQTRRVTCQKLKASGISTPVSNDMCTQVAKRPVDTQACNQQLCWEWAFSSWGQCNGPCIGPHLAVQHRQVFCQ TRDGITLPSEQSALPRPVSTQNCWSEACSVHWRVSLWTLCTATCGNYGFQSRRVEC VHARTNKAVPEHLCSWGPRPAN WQRCNITPCENMECRDTTRYCEKVKQLKLCQLSQFKSRCCGTCGKA

靶点信息

研究背景	This gene encodes a secreted protein and member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motif) family. This protein lacks the metalloproteinase and disintegrin-like domains, which are typical of the ADAMTS family, but contains other ADAMTS domains, including the thrombospondin type 1 motif. This protein may have important functions in the extracellular matrix. Alternative splicing results in multiple transcript variants encoding distinct proteins.
------	--

基因ID	92949
基因名	ADAMTSL1
Swiss	Q8N6G6
别名	ADAMTSL1;ADAMTSL-1;ADAMTSR1;C9orf94;PUNCTIN

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn