

LHX4 Rabbit pAb

货号: B16097

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	43kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 121-390 of human LHX4 (NP_203129.1).
序列	NRQLATGDEFYLMEDGRLVCKEDYETAKQNDDSEAGAKRPRTTITAKQLET LKNAYKNSPKPARHVVREQLSSETGLDMRV VQVWFQNRRAKEKRLKKDAGRHRWGQFYKSVKRSRGSSKQEKESSAEDCGVSDSELSFREDQILSELGHTNRIYGNVG DVTGGQLMNGSFMSMDGTGQSYQDLRDGSPYGIPQSPSSISLPSHAPLLNGLDYTEVDSNLGIIAHAGQGVSTLRAMAG GPTSDISTGSSVGYPDFPTSPGSWLDEMDHPPF

靶点信息

研究背景	This gene encodes a member of a large protein family which contains the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein is a transcription factor involved in the control of differentiation and development of the pituitary gland. Mutations in this gene cause combined pituitary hormone deficiency 4.
------	--

基因ID	89884
基因名	LHX4
Swiss	Q969G2
别名	LHX4;CPHD4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn