

MTHFD1L Rabbit pAb

货号: B16137

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	29kDa/105kDa
实测分子量	106kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,placenta,Mouse brain,Mouse heart,Rat brain
细胞定位	Mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 759-978 of human MTHFD1L (NP_056255.2).
序列	GVPLKKEYTEENIQLVADGCCNLQKQIQITQLFGVPVVVALNVFKTDTRAIEDLVCELAKRAGAFDAVPCYHWSVGGKGSV DLARAVREAASKRSRFQFLYDVQVPIVDKIRTIAQAVYGAKDIELSPEAQAKIDRYTQQGFGNLPICMAKTHLSLSHQPKK GVPRDFILPISDVRASIGAGFIYPLVGTMTMPGLPTRPCFYDIDLDTETEQVKGLF

靶点信息

研究背景	The protein encoded by this gene is involved in the synthesis of tetrahydrofolate (THF) in the mitochondrion. THF is important in the de novo synthesis of purines and thymidylate and in the regeneration of methionine from homocysteine. Several transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	25902
基因名	MTHFD1L
Swiss	Q6UB35
别名	MTHFD1L;FTHFSDC1;MTC1THFS;dj292B18.2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn