

ASCC3 Rabbit pAb

货号: B16141

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:100
理论分子量	13kDa/83kDa/251kDa
实测分子量	13kDa/83kDa/251kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Jurkat,HeLa,NCI-H460,Mouse liver,Mouse spleen,Mouse brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human ASCC3 (NP_071374.1).
序列	MALPRLTGALRSFSNVTKQDNYNEEVADLKIKRSLHEQVLDLGLTWKKIIFLNEKLEKSKMQSINEDLKDILHAAKQIEVN CPFQKRRLDGKEEDEKMSRASDRFRGLR

靶点信息

研究背景	This gene encodes a protein that belongs to a family of helicases that are involved in the ATP-dependent unwinding of nucleic acid duplexes. The encoded protein is the largest subunit of the activating signal cointegrator 1 complex that is involved in DNA repair and resistance to alkylation damage. Alternate splicing results in multiple transcript variants.
基因ID	10973

基因名	ASCC3
Swiss	Q8N3C0
别名	ASCC3;ASC1p200;HELIC1;RNAH

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn