

EYA3 Rabbit pAb

货号: B16173

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	48kDa/56kDa/57kDa/58kDa/62kDa
实测分子量	63kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,HT-29,Mouse liver,Mouse eye,Rat eye,Rat brain,Rat liver
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human EYA3 (NP_001981.2).
序列	MEEEQDLPEQPVKKAKMQESGEQTISQVSNPDVSDQKPETSSLASNLPMSEEIMTCTDYIPRSSNDYTSQMYSAKPYAHIL SVPVSETAYPGQTQYQTLQQTQPYAVYPQATQTYGLPPFGALWPGMKPESGLIQTPSPSQHSVLTCTTGLTTSQPSPAHY SYPIQASSTNASLISTSSTIANIPAAAVASISNQDYPTYTILGQNQYQACYPSSSFGVTGQTNSDAESTTLAATTYQSEKPSV MAPAPAAQRLSSGDPSTSPSLSQTTTPSKDIDDQSRKNMTSKNRGKRKADATSSQDS

靶点信息

研究背景	This gene encodes a member of the eyes absent (EYA) family of proteins. The encoded protein may act as a transcriptional activator and have a role during development. It can act as a mediator of chemoresistance and cell survival in Ewing sarcoma cells, where this gene is up-regulated via a micro-RNA that binds to the 3' UTR of the transcript. A similar protein in mice acts as a transcriptional activator. Alternative splicing of this gene results in multiple transcript variants.
------	--

基因ID	2140
基因名	EYA3
Swiss	Q99504
别名	EYA3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn