

ELMOD2 Rabbit pAb

货号: B16187

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	34kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,LO2
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of mouse ELMOD 2 (Q8BGF6).
序列	MFVSLWEFFYGHFFRFWMKWLLRQMTGKCELQRIFDTYGGAQRITYRIENSLTYSKNKVLQNA TRVAQSELDRCIANIMKE KNICSEKDTSFQICMRTCLLQITGYKQLYHDVENVRKKPYDSANAQHEKM

靶点信息

研究背景	This gene encodes one of six engulfment and motility (ELMO) domain-containing proteins. This gene is thought to play a role in antiviral responses. Mutations in this gene may be involved in the cause of familial idiopathic pulmonary fibrosis.
基因ID	255520

基因名	ELMOD2
Swiss	Q8IZ81
别名	ELMOD2;9830169G11Rik

产品验证



Western blot analysis of ELMOD2 expressed in U-87MG,LO2 using ELMOD2 Rabbit pAb at 1:1000. Second ary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer : 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn