

STARD3 Rabbit pAb

货号: B16209

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	48kDa/49kDa/50kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,SKOV3,MCF7,A-549,Mouse lung,Mouse brain,Mouse heart,Rat brain
细胞定位	Late endosome membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 170-445 of human STAR D3 (NP_006795.3).
序列	DFKVLPQEAEERWYLAAQVAVARGPLLFSGALSEGQFYSPPEFAGSDNESDEEVAGKKSFSAQEREYIRQGKEATAVV DQILAQEENWKFEKNNEYGDTVYTIEVPFHGKTFILKTFLPCPAELVYQEVILQPERMVLWNKTVTACQILQRVEDNTLISYD VSAGAAGGVVSPRDFVNVRRIERRRDRYLSSGIATSHSAKPPTHKYVRGENGPGGFIVLKSASNPRVCTFVWILNTDLKGR LPRYLIHQSLAATMFEFAFHLRQRISLGARA

靶点信息

研究背景	This gene encodes a member of a subfamily of lipid trafficking proteins that are characterized by a C-terminal steroidogenic acute regulatory domain and an N-terminal metastatic lymph node 64 domain. The encoded protein localizes to the membranes of late endosomes and may be involved in exporting cholesterol. Alternative splicing results in multiple transcript variants.
------	--

基因ID	10948
基因名	STARD3
Swiss	Q14849
别名	STARD3;CAB1;MLN64;es64

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn