

AHCYL1 Rabbit pAb

货号: B16211

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	53kDa/58kDa
实测分子量	59kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,U-251MG,SKOV3,Mouse brain,Rat spinal cord
细胞定位	Apical cell membrane,Cytoplasm,Endoplasmic reticulum,Microsome,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 31-260 of human AHCYL 1 (NP_006612.2).
序列	MATVTKAPKKQIQFADDMQEFTKFPTKTGRRSLRSISQSSTDSYSSAASYTDSSDDEVSPREKQQTNSKGSSNFCVKNI KQAEFGRREIEIAEQDMSALISLRKRAQGEKPLAGAKIVGCTHITAQTAVLIETLCALGAQCRWSACNIYSTQNEVAAALAEA GVAVFAWKGESEDDFWWCIDRCVNMDBGWQANMILDDGGDLTHWVYKKYPNVFKKIRGIVEESVTGVH

靶点信息

研究背景	The protein encoded by this gene interacts with inositol 1,4,5-trisphosphate receptor, type 1 and may be involved in the conversion of S-adenosyl-L-homocysteine to L-homocysteine and adenosine. Several transcript variants encoding two different isoforms have been found for this gene.
基因ID	10768

基因名	AHCYL1
Swiss	O43865
别名	AHCYL1;DCAL;IRBIT;PPP1R78;PRO0233;XPVKONA

产品验证



Western blot analysis of AHCYL1 expressed in Mouse brain,Rat spinal cord,293T using AHCYL1 Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug p er lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s

实验步骤

访问官网浏览详情: www.ablybio.cn