

ARPC3 Rabbit pAb

货号: B16214

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	20kDa
实测分子量	21kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,SKOV3,HeLa,SH-SY5Y,Mouse spleen,Mouse lung,Rat ovary
细胞定位	Cell projection,Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human ARPC3 (NP_001265485.1).
序列	MPAYHSSLMDPDTKLIGNMALLPIRSQFKGPAPRETKDTDIVDEAIYFKANVFFKNYEIKNEADRTLIIYITLYISECLKKLQKC NSKSQGEKEMYTLGITNFPPIGEPGFPLNAIYAKPANKQEDEVMRAYLQQLRQETGLRLCEKVFDPPQNDKPSKWWTCFVK RQFMNKSLSGPGQ

靶点信息

研究背景	This gene encodes one of seven subunits of the human Arp2/3 protein complex. The Arp2/3 protein complex has been conserved through evolution and is implicated in the control of actin polymerization in cells. Alternatively spliced transcript variants have been found for this gene.
基因ID	10094

基因名	ARPC3
Swiss	O15145
别名	ARPC3;ARC21;p21-Arc

产品验证



Western blot analysis of ARPC3 expressed in Mouse spleen,Rat ovary,U-251MG using ARPC3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn