

TTF2 Rabbit pAb

货号: B16219

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	75kDa/129kDa
实测分子量	130kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,HL-60,SKOV3,Jurkat
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human TTF2 (N P_003585.3).
序列	MEEVRCPEHGTFCFLKTGVRDGPNGKGSFYVCRADTCSFVRATDIPVSHCLLHEDFVVELQGLLLPQDKKEYRLFFRCIRS KAEGKRWCGSIPWQDPDSKEHSVSNKSQHASETFHHSSNWLRNPFKVLDPALWKQLIKGEGEEKKADKKQREK GDQLFDQKKEQKPEMMEKDLSSGLVPKKKQSVVQEKQEEGAEIQCEAETGGTHKRDFSEIKSQQCQGNELTRPSASS QEKSSGKSQDVQRESEPLREKVTQLLPQNVHSHN

靶点信息

研究背景	This gene encodes a member of the SWI2/SNF2 family of proteins, which play a critical role in altering protein-DNA interactions. The encoded protein has been shown to have dsDNA-dependent ATPase activity and RNA polymerase II termination activity. This protein interacts with cell division cycle 5-like, associates with human splicing complexes, and plays a role in pre-mRNA splicing.
------	--

基因ID	8458
基因名	TTF2
Swiss	Q9UNY4
别名	TTF2;HuF2;ZGRF6

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn