

TRPS1 Rabbit pAb

货号: B16221

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:2000 - 1:10000
理论分子量	141kDa/143kDa
实测分子量	160KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse kidney,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 20-260 of human TRPS1 (NP_054831.2).
序列	PPLRNVASEGEGQILEPIGTESKVGSKNKEFSADQMSENTDQSDAAELNHKEEHSLVQDPSSSSKKDLKSAVLSEKAGF NYESPSKGGNFPSFPHDEVTDNRNMLAFSSPAAGGVCEPLKSPQRAEADDPQDMACTPSGDSLETKEDQKMSPKATEET GQAQSGQANCQGLSPVSVASKNPQVPSDGGVRLNKSKTDLLVNDNPDPAPLSPELQDFKCNICGYGYGNDPTDLIKHF RKYH

靶点信息

研究背景	This gene encodes a transcription factor that represses GATA-regulated genes and binds to a dynein light chain protein. Binding of the encoded protein to the dynein light chain protein affects binding to GATA consensus sequences and suppresses its transcriptional activity. Defects in this gene are a cause of tricho-rhino-phalangeal syndrome (TRPS) types I-III.
------	--

基因ID	7227
基因名	TRPS1
Swiss	Q9UHF7
别名	TRPS1;GC79;LGCR

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn