

SLC4A2 Rabbit pAb

货号: B16226

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	135kDa/136kDa/137kDa
实测分子量	140kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,Jurkat,HT-29
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human SLC4 A2 (NP_001186623.1).
序列	DDGGASGRPLPKAQPGRHSYNLQERRRIGSMTGAEQALLPRVPTDEIEAQLATADLDLMKSHRFEDVPGVRRHLVRKN AKGSTQSGREGREPGTPRARP

靶点信息

研究背景	This gene encodes a member of the anion exchanger family of membrane transport proteins. The encoded protein regulates intracellular pH, biliary bicarbonate secretion, and chloride uptake. Reduced expression of this gene may be associated with primary biliary cirrhosis (PBC) in human patients, while differential expression of this gene may be associated with malignant hepatocellular carcinoma, colon and gastric cancers.
基因ID	6522

基因名	SLC4A2
Swiss	P04920
别名	SLC4A2;AE2;BND3L;EPB3L1;HKB3;NBND3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn