

FGF6 Rabbit pAb

货号: B16246

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	22kDa
实测分子量	22kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat skeletal muscle
细胞定位	Secreted,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FGF6 (NP_066276.2).
序列	MALGQKLFITMSRGAGRLQGTLWALVFLGILVGMVVPSPAGTRANNTLLDSRGWGTLLSRSRAGLAGEIAGVNWESGYLVGIKRQRRLYCNVGIGFHLQV

靶点信息

研究背景	The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This gene displayed oncogenic transforming activity when transfected into mammalian cells. The mouse homolog of this gene exhibits a restricted expression profile predominantly in the myogenic lineage, which suggested a role in muscle regeneration or differentiation.
------	---

基因ID	2251
基因名	FGF6
Swiss	P10767
别名	FGF6;HBGF-6;HST2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn