

CHI3L2 Rabbit pAb

货号: B16254

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	34kDa/42kDa/43kDa
实测分子量	40KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 27-390 of human CHI3L2 (NP_003991.2).
序列	YKLVCYFTNWSQDRQEPGKFTPENIDPFLCSHLIYSFASIENNKVIKDKSEVMLYQTINSLKTKNPKLKILLSIGGYLFGSKGF HPMVDSSSTRLEFINSIILFLRNHNFDGLDVSWIYPDQKENTHFTVLIHELAEAFQKDFTKSTKERLLLTAGVSAGRQMIDN SYQVEKLAKDLDFINLLSFDHGSWEKPLITGHNSPLSGWQDRGPSSYYNVEYAVGYWIIHKGMPSEKVVMGIPTYGHSE TLASAETTVGAPASGPGAAGPITESSGFLAYYEICQFLKGAKITRLQDQQVPYAVKGNQWVGYYDDVKSMETKVQFLKNLN LGGAMIWSIDMDDFTGKSCNQGPYPLVQAVKRSLGSL

靶点信息

研究背景	The protein encoded by this gene is similar to bacterial chitinases but lacks chitinase activity. The encoded protein is secreted and is involved in cartilage biogenesis. Several transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	1117
基因名	CHI3L2
Swiss	Q15782
别名	CHI3L2;CHIL2;YKL-39;YKL39

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn