

TSSK3 Rabbit pAb

货号: B16271

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	30kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Mouse brain,Mouse kidney,Rat testis
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-268 of human TSSK3 (NP_443073.1).
序列	MEDFLLSNGYQLGKTIGEGTYSKVKEAFSKKHQRKVAIKVIDKMGGPEEFIQRFLPRELQIVRTL DHKNIIQVYEMLESADG KICLVMELAEGGDVFD CVLNGGPLPESRAKALFRQMVEAIRYCHGCGVAHRDLKCENALLQGFNLKLTDFGFAKVLPKSH RELSQTFCGSTAYAAPEVLQGIPHDSKKGDVWSMGVVLVVMLCASLPFDDTDIPKMLWQQQKGVSFPTHLSISADCQDL LKRLLEPDMILRPSIEEVSWHPWLAST

靶点信息

研究背景	This gene encodes a kinase expressed exclusively in the testis that is thought to play a role in either ger m cell differentiation or mature sperm function.
基因ID	81629

基因名	TSSK3
Swiss	Q96PN8
别名	TSSK3;SPOGA3;STK22C;STK22D;TSK3

产品验证



Western blot analysis of TSSK3 expressed in Mouse testis,Rat testis,Mouse brain using TSSK3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn