

LSM11 Rabbit pAb

货号: B16302

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	39kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HepG2,MCF7,HL-60,Rat brain,Rat pancreas
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 131-360 of human LSM1 1 (NP_775762.1).
序列	RRGPGRSRKAPRNVLTRMPLHEGSPLGELHRCIREGVKVNVIHRTFKGLRGVCTGFLVAFDKFWMALTDVDETYRKPVLGKAYERDSSLTLTRLFDRLLQDSSKKEADSKSAVEDSTLSRYSTWKLASVWGRADTGRGSHKRSRVPSSLQASAREESRSELSGRTTRTDGSSVGGTFSRATTLSRGQSRKKKRKPKVDYQQVFTRHINQIFIRGENVLLVHLAQ

靶点信息

研究背景	Component of the U7 snRNP complex that is involved in the histone 3'-end pre-mRNA processing (By similarity. Increases U7 snRNA levels but not histone 3'-end pre-mRNA processing activity, when overexpressed. Required for cell cycle progression from G1 to S phases. Binds specifically to the Sm-binding site of U 7 snRNA.
基因ID	134353

基因名	LSM11
Swiss	P83369
别名	LSM11

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn