

Smad9 Rabbit pAb

货号: B16304

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC IF/ICC
推荐浓度	IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	48kDa/52kDa
实测分子量	45-52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293T,A-549
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Smad9 (NP_005896.1).
序列	MHSTTPISLSFTSPAVKRLLGWKQGDEEEKWAEKAVDSL VKLKKKKGAMDELERALSCPGQPSKCVTIPRSLDGRLQ VSHRKGLPHVIYCRVWRWPD LQSHHELKPLECCEFPFGSKQKEVCINPYHYRRVETPVLPPVLVPRHSEYNPQLSLLAKFRS ASLHSEPLMPHNATYPDSFQQPPCSALPPSPSHAFSQSPCTASYPHSPGSPSEPESPYQHSDFRPVCYEEPQHWCSVAYY ELNNRVGETFQASSRSVLIDGFTDPSNNRNRFLGLLSNVNRNSTIENTRRHIGKGVH

靶点信息

研究背景	The protein encoded by this gene is a member of the SMAD family, which transduces signals from TGF-beta family members. The encoded protein is activated by bone morphogenetic proteins and interacts with SMAD4. Two transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	4093
基因名	SMAD9
Swiss	O15198
别名	SMAD9;MADH6;MADH9;PPH2;SMAD8;SMAD8/9;SMAD8A;SMAD8B

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn