

UBA6 Rabbit pAb

货号: B16305

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:50 - 1:100
理论分子量	37kDa/43kDa/65kDa/117kDa
实测分子量	59kDa, 118kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,SKOV3,293T,Mouse brain
细胞定位	cytoplasm,cytosol,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human UBA6 (NP_060697.4).
序列	MEGSEPVAAHQGEEASCSSWGTGSTNKNLPIMSTASVEIDDALYSRQRYVLGDTAMQKMAKSHVFLSGMGGLGLEIAK NLVLAGIKAVTIHDEKCAWDLGTFNFFLSEDDVVNKRNRAEAVLKHIAELNPYVHVTSSSVPFNETTDLNFLDKYQCVVLT EMKLPLQKKINDFCRSQCPIKFISADVHGIWSRLFCDFGDEFVLDTTGEEPKEIFISNITQANPGIVTCLENHPHKLETGQF LTFREINGMTGLNGSI

靶点信息

研究背景	Modification of proteins with ubiquitin (UBB; MIM 191339) or ubiquitin-like proteins controls many signalin g networks and requires a ubiquitin-activating enzyme (E1), a ubiquitin conjugating enzyme (E2), and a u biquitin protein ligase (E3). UBE1L2 is an E1 enzyme that initiates the activation and conjugation of ubiqui tin-like proteins (Jin et al., 2007 [PubMed 17597759]).[supplied by OMIM, Mar 2008]
------	---

基因ID	55236
基因名	UBA6
Swiss	A0AVT1
别名	UBA6;E1-L2;MOP-4;UBE1L2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn