

PPP3R2 Rabbit pAb

货号: B16312

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa
实测分子量	16kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,H9C2,Mouse heart,Rat heart
细胞定位	cytosol,mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-173 of human PPP3R2 (NP_671709.1).
序列	MSTMGNEASYPAEMCSHFDNDEIKRLGRRFKKLDLDKSGSLSVEEFMSLPELRHNPLVRRVIDVFDTDGDGEVDFKEFIL GTSQFSVKGDEEQKLRFAFSIYDMDKDGYSISNGELFQVLKMMVGNNLTDWQLQQLVDKTIILDKDGDGKISFEESAVV RDLEIHKKLV LIV

靶点信息

研究背景	Regulatory subunit of calcineurin, a calcium-dependent, calmodulin stimulated protein phosphatase. Confers calcium sensitivity.
基因ID	5535

基因名	PPP3R2
Swiss	Q96LZ3
别名	PPP3R2;PPP3RL

产品验证



Western blot analysis of PPP3R2 expressed in Mouse heart,Rat heart,HeLa using PPP3R2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn