

OX40 Rabbit pAb

货号: B16316

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	29kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,BT-474,SW620,MCF7,Mouse thymus,Mouse skeletal muscle,Mouse lung,Mouse intestine,Rat liver
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 29-214 of human OX40 (NP_003318.1).
序列	LHCVGDTYPSNDRCCHECRPGNGMVSRCRSQNTVCRPCGPGFYNDVVSSKPKPCTWCNLRSGSERKQLCTATQDTV CRCRAGTQPLDSYKPGVDCAPCPGHHFSPGDNQACKPWTNCTLAGKHTLQPASNSSDAICEDRDPPATQPQETQGPPA RPITVQPTEAWPRTSQGPSTRPVEVPGGRA

靶点信息

研究背景	The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor has been shown to activate NF-kappaB through its interaction with adaptor proteins TRAF2 and TRAF5. Knockout studies in mice suggested that this receptor promotes the expression of apoptosis inhibitors BCL2 and BCL2L1/BCL2-XL, and thus suppresses apoptosis. The knockout studies also suggested the roles of this receptor in CD4+ T cell response, as well as in T cell-dependent B cell proliferation and differentiation.
------	--

基因ID	7293
基因名	TNFRSF4
Swiss	P43489
别名	TNFRSF4;ACT35;CD134;IMD16;OX40;TXGP1L

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn