

CAPN5 Rabbit pAb

货号: B16345

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	73kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,Mouse thymus,Mouse brain,Mouse spinal cord
细胞定位	cell surface,cytoplasm,extracellular exosome,focal adhesion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 361-640 of human CAPN5 (NP_004046.2).
序列	GAWTLHEDPRQNRGGGCINHKDTFFQNPQYIFEVKKPEDEVLCIQRPKRSTRREGKGENLAIGFDIYKVEENRQYRMHS LQHKAASSIYINRSVFLRTDQPEGRYVLIPTTFEPGHTGEFLLRVFTDVPSNCRELRLDEPPHTCWSSLCGYPQLVTQVHVL GAAGLKDSPTGANSYVIKCEGDKVRSQVQKGTSTPEYNVKGIFYRKKLSQPITVQVWNHRVLKDEFLGQVHLKADPDNL QALHTLHLDRNSRQPSNLPGTAVVHLSSTSLMAV

靶点信息

研究背景	Calpains are calcium-dependent cysteine proteases involved in signal transduction in a variety of cellular processes. A functional calpain protein consists of an invariant small subunit and 1 of a family of large subunits. CAPN5 is one of the large subunits. Unlike some of the calpains, CAPN5 and CAPN6 lack a calmodulin-like domain IV. Because of the significant similarity to Caenorhabditis elegans sex determination gene tra-3, CAPN5 is also called as HTRA3.
------	--

基因ID	726
基因名	CAPN5
Swiss	O15484
别名	CAPN5;ADNIV;HTRA3;VRNI;nCL-3;calpain-5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn