

CDYL2 Rabbit pAb

货号: B16355

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	56kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,U-937,HeLa,Mouse brain,Mouse spleen,Mouse testis,Mouse lung,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 60-230 of human CDYL2 (NP_689555.2).
序列	MSKDKRIKSGKQSSTSKLLRDSRGPSVEKLSHRPSDPGKSKGTSHKRKRINPPLAKPKKGYSGKPSSGGDRATKTVSYRT TPSGLQIMPLKKSQNGMENGDAGSEKDERHFGNGSHQPGLDLNDHVGEQDMGECDVNHATLAENGLGSALTNGGLN LHSPVKRKLEAEKDY

靶点信息

研究背景	Predicted to enable transcription corepressor activity. Predicted to be involved in negative regulation of tr anscription, DNA-templated. Predicted to be active in nucleus.
基因ID	124359

基因名	CDYL2
Swiss	Q8N8U2
别名	CDYL2;PCCP1

产品验证



Western blot analysis of CDYL2 expressed in Mouse brain,Rat brain,MCF7 using CDYL2 Rabbit pAb at 1: 1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn