

NAA60 Rabbit pAb

货号: B16356

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	19kDa/20kDa/27kDa/28kDa
实测分子量	27KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293F
细胞定位	Golgi apparatus membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-242 of human NAA60 (NP_079121.1).
序列	MTEVVPSSALSEVSLRLLCHDDIDTVKHLCGDWFPiEYPDSWYRDITSNKKFFSLAATYRGAIVGMIVAEIKNRTKIHKEDG DILASNFSVDTQVAYILSLGVVKEFRKHGIGSLLLLESKDKHISTTAQDHCKAIYLVLTNTNTAINFYENRDFKQHHYLPYYYS IRGVLKDGFTYVLYINGGHPPWTILDYIQHLGSALASLSPCSIPHRVYRQAHSLLCSFLPWSGISSKSGIEYSRTM

靶点信息

研究背景	This gene encodes an enzyme that localizes to the Golgi apparatus, where it transfers an acetyl group to the N-terminus of free proteins. This enzyme acts on histones, and its activity is important for chromatin assembly and chromosome integrity. Alternative splicing and the use of alternative promoters results in multiple transcript variants. The upstream promoter is located in a differentially methylated region (DMR) and undergoes imprinting; transcript variants originating from this position are expressed from the maternal allele.
基因ID	79903
基因名	NAA60
Swiss	Q9H7X0
别名	NAA60;HAT4;NAT15;NatF;hNaa60

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn