

# TNFRSF21 Rabbit pAb

货号: B16367

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                        |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB                                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                              |
| 理论分子量 | 71kDa                                                                                                  |
| 实测分子量 | 110kDa                                                                                                 |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | HT-29,A-549,SW480,SW620,Mouse kidney,Mouse brain,Mouse lung,Rat brain,Rat kidney                       |
| 细胞定位  | Cell membrane,Single-pass type I membrane protein                                                      |
| 纯化    | Affinity purification                                                                                  |

抗原信息


|      |                                                                                                                                                                                                                                                                                                            |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 371-655 of human TNFRSF21 (NP_055267.1).                                                                                                                                                                                     |
| 序列   | RKSSRTLKKGPRQDPSAIVEKAGLKKSMTPTQNREKWIYYCNGHGIDILKLVAQVGSQWKDIYQFLCNASEREVAAFSN<br>GYTADHERAYAALQHWITIRGPEASLAQLISALRQHRRNDVVEKIRGLMEDTTQLETDKLALPMSPSPSPSPSPNAKLEN<br>SALLTVEPSPQDKNKGFFVDESEPLLRC DSTSSGSSALSRNGSFITKEKKD TVLRQVRLDPCDLQPIFDDMLHFLNPEELR<br>VIEEIPQAEDKLDRLF E IIGVKSQEASQTLLDSVYSHLPDLL |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the tumor necrosis factor receptor superfamily. The encoded protein activates nuclear factor kappa-B and mitogen-activated protein kinase 8 (also called c-Jun N-terminal kinase 1), and induces cell apoptosis. Through its death domain, the encoded receptor interacts with tumor necrosis factor receptor type 1-associated death domain (TRADD) protein, which is known to mediate signal transduction of tumor necrosis factor receptors. Knockout studies in mice suggest that this gene plays a role in T-helper cell activation, and may be involved in inflammation and immune regulation. |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                           |
|-------|---------------------------|
| 基因ID  | 27242                     |
| 基因名   | TNFRSF21                  |
| Swiss | O75509                    |
| 别名    | TNFRSF21;BM-018;CD358;DR6 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)