

ZNF195 Rabbit pAb

货号: B16373

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	10kDa/64kDa/69kDa/70kDa/72kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,SW620,A-549,Mouse brain,Mouse heart,Mouse skeletal muscle,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ZNF195 (NP_001123991.1).
序列	MTLLTFRDVAIEFSLEEWKCLDLAQONLYRDVMLENYRNLFVGLTVCKPGLITCLEQRKEPWNVKRQEAADGHPEMGFH HATQACLELLGSSDLPASASQSAGITGVNHRAQPGLNVSVDKFTAMSSHFTQDLLPEQGIQDAFPKRILRGYGNCGLDNL YLRKDWESLDECKLQKDYNGLNQCSSTTHSKIFQYNKYVKIFDNFSNLHRRNISNTGEKPFKCQECGKSFQMLSFLTEHQ KIHTGKKFQK

靶点信息

研究背景	This gene encodes a protein belonging to the Krueppel C2H2-type zinc-finger protein family. These family members are transcription factors that are implicated in a variety of cellular processes. This gene is located near the centromeric border of chromosome 11p15.5, next to an imprinted domain that is associated with maternal-specific loss of heterozygosity in Wilms' tumors. Alternative splicing results in multiple transcript variants.
------	---

基因ID	7748
基因名	ZNF195
Swiss	O14628
别名	ZNF195;HRF1;ZNFP104

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn