

UBE2H Rabbit pAb

货号: B16375

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	17kDa/20kDa
实测分子量	21-42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,BT-474,293T,U-251MG,SKOV3,Mouse thymus,Mouse liver,Mouse lung,Rat lung
细胞定位	cytosol,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human UBE2H (NP_003335.1).
序列	MSSSPSPGKRRMDTDVVKLIESKHEVTILGGLNEFVVKFYGPQGTPYEGGVWKVRVDLPDKYPFKSPSIGFMNKIFHPNIDE ASGTVCLDVINQWTALYDLTNIFESFLPQLLAYPNPIDPLNGDAAAMYLHRPEEYKQKIKEYIQKYATEEALKEQEEGTGDS SSESSMSDFSEDEAQDMEL

靶点信息

研究背景	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or s hort-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-ac tivating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. Thi s gene encodes a member of the E2 ubiquitin-conjugating enzyme family. The encoded protein sequence is 100% identical to the mouse homolog and 98% identical to the frog and zebrafish homologs. Three alte rnatively spliced transcript variants have been found for this gene and they encode distinct isoforms.
------	--

基因ID	7328
基因名	UBE2H
Swiss	P62256
别名	UBE2H;E2-20K;GID3;UBC8;UBCH;UBCH2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn