

ATMIN Rabbit pAb

货号: B16406

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	72kDa/88kDa
实测分子量	68-88kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Mouse testis,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 160-360 of human ATM N (NP_056066.2).
序列	HAEEKHKCSKCSNSYGTWDLKRHAEDCGKTFRCTGCPYASRTALQSHIYRTGHEIPAHRDPPSKKRKMENCAQNQK LSNKTIESLNNQPIRPDTQELESEIKLEPSFEDSCGSNTDKQTLTTPPRYPQKLLLPKPKVALVKLPVMQFSVMPFVPTA DSSAQPVVLGVDQGSATGAVHLMPLSVGTLILGLDSEAC

靶点信息

研究背景	Transcription factor. Plays a crucial role in cell survival and RAD51 foci formation in response to methylating DNA damage. Involved in regulating the activity of ATM in the absence of DNA damage. May play a role in stabilizing ATM. Binds to the DYNLL1 promoter and activates its transcription.
基因ID	23300

基因名	ATMIN
Swiss	O43313
别名	ATMIN;ASCIZ;ZNF822

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn