

GJA5 Rabbit pAb

货号: **B16422**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	40kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse heart,Rat heart
细胞定位	Cell junction,Cell membrane,Multi-pass membrane protein,gap junction
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 229-358 of human GJA5 (NP_859054.1).
序列	LGWKKIRQRFVKPRQHMAKCQLSGPSVGIVQSCTPPPDFNQCLENGPGGKFFNPFSNNMASQQNTDNLVTEQVRGQE QTPGEGFIQVRYGQKPEVPNGVSPGHRPLPHGYHSDKRRRLSKASSKARSDDLVS

靶点信息

研究背景	This gene is a member of the connexin gene family. The encoded protein is a component of gap junctions , which are composed of arrays of intercellular channels that provide a route for the diffusion of low molecular weight materials from cell to cell. Mutations in this gene may be associated with atrial fibrillation. Alternatively spliced transcript variants encoding the same isoform have been described.
基因ID	2702

基因名	GJA5
Swiss	P36382
别名	GJA5;ATFB11;CX40

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn