

RSRC1 Rabbit pAb

货号: **B16431**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	31kDa/38kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,Mouse brain
细胞定位	Cytoplasm,Nucleus,Nucleus speckle
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 117-276 of human RSRC1 (NP_001258763.1).
序列	EHLPPAEQAKARLQLVLEAAAKADEALKAKERNEEEAKRRKEEDQATLVEQVKRVKEIEAIESDSFVQQQTFRSSKEVKKSVE PSEVKQATSTSGPASAVADPPSTEKEIDPTSIPTAIKYQDDNSLAHPNLFIEKADAEEKWFKRLIALRQERLMGSPVA

靶点信息

研究背景	This gene encodes a member of the serine and arginine rich-related protein family. The encoded protein is involved in both constitutive and alternative mRNA splicing. This gene may be associated with schizophrenia. A pseudogene of this gene is located on chromosome 9. Alternative splicing results in multiple transcript variants encoding different isoforms.
------	--

基因ID	51319
基因名	RSRC1
Swiss	Q96IZ7
别名	RSRC1;BM-011;SFRS21;SRrp53

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn