

KIAA1456 Rabbit pAb

货号: B16438

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	14kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,BT-474,Mouse pancreas
细胞定位	cytoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-330 of human KIAA 1456 (NP_065895.2).
序列	VIHHFSTKQRRIRAIKEMARVLVPGGQLMIYVWAMEQKNRHFQKQDVLVPWNRALCSQLFSESSQSGRKRQCGYPERGH PYHPPCSECSVCFKQCGSKRSHSVGYEPAMARTCFANISKEGEEEEYGFYSTLGKSFRSWFFSRSLDESTLRKQIERVR PLKNTEVWASSTVTVQPSRHSSLDQEPFSTKGQSLDEEVFVESSSGKHLEWLRAPGTL

靶点信息

研究背景	May modifie wobble uridines in specific arginine and glutamic acid tRNAs. Acts as a tumor suppressor by promoting the expression of LIN9.
基因ID	57604

基因名	KIAA1456
Swiss	Q9P272
别名	KIAA1456;C8orf79;TRM9L

产品验证



Western blot analysis of KIAA1456 expressed in HepG2,BT-474,Mouse pancreas using KIAA1456 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: www.ablybio.cn