

PLA2G2D Rabbit pAb

货号: B16448

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:10 - 1:100
理论分子量	16kDa
实测分子量	22KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Daudi,RPMI 8226
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 21-145 of human PLA2G 2D (NP_036532.1).
序列	GILNLNKMVKQVTGKMPILSYWPYGCHCGLGGRGQPKDATDWCCQTHDCCYDHLKTQGCSIYKDYYRYNFSQGNIHCS DKGSWCEQQLCACDKEVAFCLKRNLDTYQKRLRFYWRPHCRGQTPGC

靶点信息

研究背景	This gene encodes a secreted member of the phospholipase A2 family, and is found in a cluster of related family members on chromosome 1. Phospholipase A2 family members hydrolyze the sn-2 fatty acid ester bond of glycerophospholipids to produce lysophospholipids and free fatty acid. This gene may be involve d in inflammation and immune response, and in weight loss associated with chronic obstructive pulmonar y disease. Alternative splicing results in multiple transcript variants encoding different isoforms.
------	--

基因ID	26279
基因名	PLA2G2D
Swiss	Q9UNK4
别名	PLA2G2D;PLA2IID;SPLASH;sPLA2-IID;sPLA2S

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn