

TRMT1 Rabbit pAb

货号: B16464

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	69kDa/72kDa
实测分子量	90kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,HL-60,Mouse heart
细胞定位	nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 460-659 of human TRM T1 (NP_060192.1).
序列	PSLLQLRSALLHADFRVSLSHACKNAVKTDAPASALWDIMRCWEKECPVKRERLSETSPAFRILSVEPRLQANFTIREDANPSSRQRGLKRFQANPEANWGP RP RP RP PGGKADEAMEERRRLLQNK RKEPPEDVAQRAARLKTFPCKRFKEGTCQRGDQCCYSHSPPTPRVSADAAPDCPETSNTPPGPGAAAGPGID

靶点信息

研究背景	Dimethylates a single guanine residue at position 26 of most tRNAs using S-adenosyl-L-methionine as donor of the methyl groups.
基因ID	55621

基因名	TRMT1
Swiss	Q9NXH9
别名	TRMT1;TRM1

产品验证



Western blot analysis of TRMT1 expressed in MCF7,HL-60,Mouse heart using TRMT1 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Bl ocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn